



ARTISAN[®] Aphakia


#1 Backup Lens in complicated cataract cases

Based on the long term experience of Iris Fixation, the ARTISAN[®] Aphakia IOL is a predictable, safe, high precision implant, that corrects the eye when it is not correctable by other means.

Pediatric lenses are available upon request for juvenile or congenital cataract cases.



// Specifications Model 205



Optic:	5.0 mm, biconvex*
Body:	5.4 mm
Overall diameter:	8.5 mm
A-Constant	115.0 (A Scan) 115.7 (IOL Master; SRK T) 115.7 (IOL Master; SRK II) -0.160 (IOL Master; Haigis a0) 0.400 (IOL Master; Haigis a1) 0.100 (IOL Master; Haigis a2) 3.62 (IOL Master; Hoffer-Q pACD) -0.08 (IOL Master; Holladay 1 sf)
Dioptric Powers:	2.0 D to 30.0 D (1.0 D increments) 14.5 D to 24.5 D (0.5 D increments)
Also available:	ARTISAN [®] Pediatric Aphakia 4.4/6.5 and 4.4/7.5 Designed for small eyes

*+2.0 D to +9.0 D Convex-plano

// Main Features & Benefits

- ✓ Iris fixation
- ✓ One size fits all
- ✓ Long term clinical experience
- ✓ Predictable , stable, reliable
- ✓ Also fit for retro pupillary fixation

// ARTISAN[®] Training Course

To obtain the best results, OPHTEC continuously conducts specialized training courses throughout the world to develop and maintain surgical skills. OPHTEC's condition to use the product correctly is to provide ARTISAN[®] IOLs only to certified professionals. For more information on training courses send an e-mail to marketing@ophtec.com, or call +31 50 525 1944. Go to www.ophtec.com for upcoming training courses.

// Toric ARTISAN[®] Aphakia IOL

The Toric Artisan Aphakia IOL is a Custom Made IOL. Contact your representative for specific powers and more information.

