



ArtiFix Holding Forceps

for all ARTISAN® Family IOLs



Features & Benefits

- Reduces the risk of losing the depth of the anterior chamber
- Holds all Artisan® lenses very stable: No rotation of the lens
- Minimal enlarging of the incision
- Fit for standard implantation and retro pupillary fixation of the Artisan® Aphakia



Stable hold; no rotation of the lens.



Correct holding position for all ARTISAN® lenses; minimal enlarging of the incision.

ARTISAN® Aphakia IOL	ARTISAN® Phakic IOLs
• Standard implantation* • Retro pupillary fixation*	• ARTISAN® Myopia • ARTISAN® Hyperopia • ARTISAN® Toric

* See Artisan Aphakia IOL Training Manual for details of the implantation technique

